



 **APELEM**

PRÉSENTE

CHEST & BONE ROOMS




CAMARGUE
SERIES

 **DMS**
IMAGING
A DIVISION OF THE DMS GROUP

OUR CHEST & BONE ROOMS

CAMARGUE series: Full range of radiographic systems with high performances for all purpose of chest and bone exams



**CHEST & BONE
ROOM**




Compatible with:

- ✓ All cassette-size detectors up to 43 x 43 cm
- ✓ Fixed flat panel detector (cassette-size 43 x 43 cm)

OUR CHEST & BONE ROOMS


CAMARGUE series: 5 configurations to meet all needs



**CHEST & BONE
ROOM**



- **BASIC SYSTEM – MANUAL**



- **MOTORIZED SYSTEMS – FLOOR TUBE STAND**
AUTO-TRACKING, AUTO-STITCHING



- **MOTORIZED SYSTEMS – CEILING TUBE STAND**
AUTO-TRACKING, AUTO-STITCHING



BASIC SYSTEM - MANUAL

CAMARGUE

MODEL: Z7J32003

CAMARGUE

MODEL: Z7J32003



- **Basic system, manual**
- **Reduced footprint**
- **Floor or Ceiling Tube Stand**
- **Table with fixed height, elevating (optional)**
- **Non-tilting wall stand**
- **Oscillating and non-removable grids**
- **Optional RT movements**
- **No tomography**



MOTORIZED SYSTEMS – FLOOR TUBE STAND

CAMARGUE

2 MODELS: Z7J32004 & Z7J32005

CAMARGUE

MODEL: Z7J32004



- **Motorized system - Auto-tracking**
- **Floor Tube Stand**
- **Tube stand: vertical motorized movement**
- **Elevating table (45-90 cm)**
- **Wall stand – optional tilting**
- **High-resolution, fixed et removable grids**
- **Membrane with digital display**
- **Vertical Auto-tracking with table and wall stand**

CAMARGUE

MODEL: Z7J32004

- **Vertical Auto-tracking with table and wall stand**

- ❑ Tube: Follower
- ❑ Table and wall stand: References



- **Horizontal Auto-tracking (auto-centering) - Optional**

- ❑ Tube: Reference
- ❑ Bucky in the table: Follower



CAMARGUE

MODEL: Z7J32005




- **Motorized system - Auto-stitching**
- **Floor Tube Stand**
- **Elevating table (45-90 cm)**
- **Wall stand – optional tilting**
- **High-resolution, fixed et removable grids**
- **10" Touchscreen**
- **Vertical Auto-tracking with table and wall stand**
- **Vertical Auto-stitching**

CAMARGUE

MODEL: Z7J32005

- **Auto-stitching vertical:** Motorized tube rotation around horizontal axis, Wall stand motorized vertical movement, stitching stand





MOTORIZED SYSTEMS – CEILING TUBE STAND

CAMARGUE

2 MODELS: Z7J32006 & Z7J32007

CAMARGUE

MODEL: Z7J32006



- **Motorized system - Auto-tracking**
- **Ceiling Tube Stand**
- **Tube stand: vertical motorized movement**
- **Elevating table (45-90 cm)**
- **Wall stand – optional tilting**
- **High-resolution, fixed et removable grids**
- **Membrane with digital display**
- **Vertical Auto-tracking with table and wall stand**

CAMARGUE


MODEL: Z7J32007



- **Motorized system - Auto-stitching**
- **Ceiling Tube Stand**
- **Elevating table (45-90 cm)**
- **Wall stand – optional tilting**
- **High-resolution, fixed et removable grids**
- **10" Touchscreen**
- **Vertical Auto-tracking with table and wall stand**
- **Vertical Auto-stitching** (Motorized tube rotation around horizontal axis, Wall stand motorized vertical movement, stitching stand)

OUR CHEST & BONE ROOMS

Examples of accessories



**CHEST & BONE
ROOM**



Wall stand patient
stretch grip



Paper roll feeder



Table lateral cassette
holder



Remote controller



ADDRESS

9 avenue du canal Philippe Lamour
30660 Gallargues-le-Montueux
— FRANCE

T > +33 4 67 50 49 00

F > +33 4 67 50 49 09