

PM-3GYN GYNECOLOGICAL EXAMINATION TABLE WITH 3 MOTORS



KEY FEATURES

- Designed for gynecological examination and operation.
- Body and Base sections are isolated against any type of fluid.
- Base section is protected against corrosion with electrostatic paint.
- Antibacterial, water-proof, fabric based surface cover, customizable color.
- 3 Motors allow the motions smoothly in a vibration-free manner:
 - Angular position of the backrest
 - Adjustable table height
 - Trendelenburg / Reverse trendelenburg
- Removable footrest.
- 4 height-adjustable feet.
- Leakage current relay protects the patient from potential electric shock
- 2 years warranty.

STANDARD ACCESSORIES

- 2 Armrests
- 2 Polyurethane back-knee rests.
- 1 Stainless steel waste bin
- 1 IV Pole

OPTIONAL ACCESSORIES

- Doctor stool
- Examination light
- Anesthesia frame



TECHNICAL SPECIFICATIONS

Height adjustment range	65 – 90 cm
Backrest adjustment range (wrt seating section)	0 – 60°
Seating trendelenburg range (wrt horizontal axis)	(-18) – 0°
Footrest adjustable angles	(-90) & 0°
Operating voltage	~AC 220-240V 50 Hz
Soft area dimensions in bed configuration (LxW)	1620 x 600 mm
Weight capacity	250 kg